



# Ovarian Shield Set

684120



- Set of the three (Infant, Youth and Adult)
- Designed especially for ovarian procedures
- 1.0mm LE
- Easy to use and explain to your patient
- Available in our Classic Nylon Colors
- Durable
- Easily wipes clean

Model	Description
684120	Complete Set of 3 female ovarian shields 3" x 3½" Infant Size Ovarian Shields 3 ¾" x 4 ¾" Youth Size Ovarian Shields 4 ¾" x 5 ½" Adult Size Ovarian Shields



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

**INFAB**