



# Triangular Gonad Shields

GND



- Set of three or available individually: Infant, Youth and Adult
- 1/2" tie strap
- Easy fastening around the waist
- .50mm LE
- Available in our Classic Nylon Colors

Model	Description
GND-S1	4" x 5.5" Infant Size Gonad Shield
GND-M1	7" x 8" Youth Size Gonad Shield
GND-L1	8" x 10" Adult Size Gonad Shield
GND-X1	Set of Gonad Shield (Infant, Youth and Adult sizes included)



REV002

1040 Avenida Acaso  
Camarillo, CA 93012, USA

tel: 805.987.5255  
fax: 805.482.8424

info@infabcorp.com  
www.infabcorp.com

**INFAB**