



# Mini C-Arm Drape

## 4003



Fits mini-size C Arms, Xi-Scan, Fluoroscanner 3" or 4" FOV. Made from durable, high quality polyethylene materials with abrasion resistance to protect from rips or tears.

<b>Highlights:</b>	Sterile Product and Latex Free
<b>Quantity:</b>	40 Drapes (20 drapes per box, 2 boxes included)
<b>Dimensions:</b>	54" x 63"
<b>Packaging Size:</b>	10" x 17" x 24"
<b>Total Weight:</b>	18 lbs.