






ACCESSORIES

Thyroid Collar without Binding

NST



Dimensions | Hook and Loop & Magnetic Closures



- Comfortable fit
- Large protection area
- Optional monogram

Model Number	Description
--------------	-------------

NST	Thyroid Collar without Binding
-----	--------------------------------

Sizes	X-Small, Standard, Medium
-------	---------------------------

Colors	Revolution 700 series & All Classic Nylon Colors
--------	--

Core	Lightweight Lead KIARMOR Bi-Layer Lead-Free
------	--

Protection	.25mm LE, .35 mm LE, .50mm LE
------------	-------------------------------

Closure	Hook and Loop Closure, Buckle, Magnet
---------	---------------------------------------

Fabric	Revolution or Classic
--------	-----------------------

Monogram	Block or <i>Script</i>
----------	-------------------------------

Model Numbers	Description
---------------	-------------

NST-XS	Thyroid Collar without Binding, Size X-Small, Hook and Loop Closure
--------	---

NST-XS-M	Thyroid Collar without Binding, Size X-Small, Magnetic Closure
----------	--

NST-XS-B	Thyroid Collar without Binding, Size X-Small, Buckle Closure
----------	--

NST-S1	Thyroid Collar without Binding, Size Standard, Hook and Loop Closure
--------	--

NST-S1-M	Thyroid Collar without Binding, Size Standard, Magnetic Closure
----------	---

NST-S1-B	Thyroid Collar without Binding, Size Standard, Buckle Closure
----------	---

NST-M1	Thyroid Collar without Binding, Size Medium, Hook and Loop Closure
--------	--

NST-M1-M	Thyroid Collar without Binding, Size Medium, Magnetic Closure
----------	---

NST-M1-B	Thyroid Collar without Binding, Size Medium, Buckle Closure
----------	---



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel : 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

