# 2023 MINI CATALOG











## PEACE OF MIND

## **In Radiation Protection**

INFAB's R&D team continually strives to safeguard healthcare professionals from ionizing scatter radiation through providing top-of-the-line protective gear. Our commitment is to give healthcare professionals the assurance and peace of mind that comes from using our reliable and expertly crafted products.



## **Empowering Customers** with Knowledge

INFAB is dedicated to providing our customers with accurate and comprehensive information about the materials and technology used in their radiation protective gear. We believe that this knowledge is crucial for making informed purchasing decisions and for understanding the level of protection being used.











## **Expert Support for Radiation Protection**

Our dedicated global team of radiation protection experts is ready to provide friendly and efficient support throughout your experience with us. With years of industry experience, our specialists are well-equipped to assist with any questions you may have.



## APRONS

Our aprons are designed with the latest technology to protect you from radiation exposure and made from high-quality materials to ensure maximum durability and comfort.

Our Revolution line offers you a variety of styles and protection levels. You can add logos, monograms and design patterns to make it your own and reflect your unique style.

### **PARTIAL OVERLAP:**

Full protection in each panel and double protection where the panels overlap in the front.

### **COLOR BLOCK DESIGN:**

Exclusive Revolution color combination. Customization options available.

### **REVOLUTION FABRIC:**

Lightweight, waterproof, abrasion resistant, flame resistant, and bacterial resistant.

### **SMART TRACK:**

Web-based garment inspection tracking system.





ghtweight Lead. KIARMOR is the only lead-free and system and sees all three universally recognized testing procedures: ASTM 1094-14, IEC 61331-1:2014, DIN 6857-1.





L103 Revolution Lumbar Vest & Skirt 711 White







303 Revolution **Velcro Front** 814 Jewelry Box Blue/ White



403 Revolution Elastic **Back Front** 615 Black/Orange



503 Revolution Front Black Belt 504 Pink/Black







703 Revolution Full Overlap Lumber Vest & Skirt 505 Navy/Black





## REVOLUTION BASE

## **APRONS**

Revolution Base provides the same comfort and protection as our premium Revolution line, but in classic nylon colors and patterns. It does not include Revolution Fabric, Color-Blocking, or Cool-Wear.



BL103

Revolution Lumbar
Vest & Skirt Base
SH-23 Art Deco

Revolution Vest & Skirt Base MIR280-09 Galaxy Black



B203
Revolution Full
Wrap Black Belt
Base
217 Mardi Gras

Revolution Front Base 207 Blue Flames





B403
Revolution Front
Elastic Back Base
184 Pink Camo







B603

Revolution Elastic
Back Full Wrap Base
SH-20 New Flames



**REV-TC** Visor Style **Thyroid Collar** 

Velcro or Buckle Closure

**NST Thyroid Collar** without Binding

Velcro, Buckle or Magnet Closure





SVT Visor Style Thyroid Collar

Velcro or Buckle Closure

TCV **Thyroid Collar** with Binding

Velcro, Buckle or Magnet Closure



## **THYROID** COLLARS

INFAB offers thyroid collars to protect against scatter radiation. Our collars can be attached to lead aprons. Choose from buckle, Velcro, or magnetic closures. Personalize with colors, patterns, monograms and logos.



## CLASSIC APRONS

The Infab Classic line of radiation protection lead aprons are available in a variety of styles. Choose from thousands of colors and patterns to make it your own. The aprons are made with standard lightweight nylon fabric and are fully customizable with colors, patterns, logos, and monograms.

### STYLES:

Choose from thousands of colors and patterns.

### MAKE IT YOUR OWN:

Fully customizable with colors, patterns, logos, and monograms.

### **NYLON FABRIC:**

Standard lightweight nylon fabric.

### **SMART TRACK:**

Web-based garment inspection tracking system.



All aprons are available in KIARMOR Bi-Layer Lead-Free and .ightweight Lead. KIARMOR is the only lead-free material that basses all three universally recognized testing procedures: ASTM :3094-14. IFC 61331-1:2014. DIN 6857-1.





MWA Special Wrap SH-7 Metallic Splatter

Full Coat Wrap

120 Camouflage



**WBT** 

2" Wide Belt SH-12 Unicorn



SDA Surgery Drop Off Shoulder Release

SH-10 Cats



**FCW** 



AWC Weight Reliever SH-2 Blocks



**EBA Elastic Black** MIR280-09 Galaxy Black

## GLASSES

## **Style + Protection**











GUCCI GG0343

**NIKE BRAZEN FURY** 

**NIKE PREMIER** 







**NIKE 7117** 

NIKE 7118

**NIKE MAVERICK RGE** 









**FITOVERS** 

**INCREDIBLES** 

EV0-01

**RADIATION PROTECTION** 




## Shields & Drapes for Interventional Procedures

Prolonged low-dose exposure to scatter radiation can increase cancer risk for healthcare professionals. Scatter Armor shields and drapes reduce this risk by minimizing overall exposure.

"The scattered radiation from the patient comprises the main source of radiation dose to staff."







## **Scatter Protection**

Scatter Armor Shields provide scatter radiation protection for radial and many other types of procedures. KIARMOR Bi-Layer is the most protective core material against scatter radiation.

### SA-MPS Scatter Armor Multipurpose Shield

Transradial Access: Left and Right Radial, Shunt, Biliary Drainage, Spinal, Fibroid Embolization.






### **Absorbent**

Scatter Armor Shields utilize the highest absorbency covering material in the industry, critical in heavy fluid situations.

## SA-MPF Scatter Armor Multipurpose Shield with Fenestration

Needle Biopsy or Puncture Procedures.



## SA-BSS Scatter Armor **Biliary Split Shield**

Biopsy Procedures, Transjugular Intraheptic Portosystemic Shunt (TIPS).




## **Flexible**

Scatter Armor Shields easily mold to the patients body for ease of use and maximum protection.


### SA-EPS Scatter Armor **EP Shield**

Subclavian Access: Bi-Ventricular Pacemaker Implants, ICS procedures, and Bi-Ventricular Pacing.



### SA-PER Scatter Armor Peripheral Shield

Peripheral procedures where the access is through the femoral artery.



## One-time Use


Scatter Armor Shields are lead-free and disposable.

## Sterile

Packaged individually in sterile packets.

## SA-AGS Scatter Armor **Angiography Shield**

Transfemoral Access: Angiography, Coronary Catheterization.



## SA-DDD Scatter Armor Dialysis Drape

Dialysis Access and Declotting Procedures.



## LOCKERS

## **Security + Protection**

Our Revolution Lockers combine protection and security for your garments. They include a passcode encrypted keypad and durable frame for secure storage. The battery-operated passcode system allows for up to 3 admin codes.




REV-LOCKER-20 Revolution 20 Hanger Mobile Apron Locker

## **REVOLUTION MOBILE APRON LOCKERS**

(Black Shown)

**REV-LOCKER-5** Revolution 5 Hanger Mobile Apron Locker

- Features digital passcode locks with keypad
- Set up to 3 administrative and user passcodes
- Battery operated lock with backup battery and emergency power supply



**REV-LOCKER-10** Revolution 10 Hanger Mobile Apron Locker

## BARRIERS

## **Mobility & Protection**

Protect yourself from head-to-toe with any of our durable mobile barriers. Most of our barriers come equipped with the option of either an acrylic or a lead shield. Optional color, metal type, height, width, and window configurations available. INFAB's mobile barriers are guaranteed to make your life easier with its convenient mobility, and easy-to-clean features.







## **BARRIERS**

\*All Glass comes in Lead or Acrylic (-AC)



## **683492** Large Window Mobile Barrier

SIZE: 801/4" H x 401/2" W x 25" D WINDOW: 72" H x 36" W

## **683476** X-Tall Mobile Barrier

SIZE: 78¼" H x 40½" W x 25"D WINDOW: 60" H x 36" W

## 683488 X-Wide Mobile Barrier

SIZE: 74½" H x 76½" W x 25"D WINDOW: 36" H x 72" W

## **683473** Mid-Tall Mobile Barrier

SIZE: 74½" H x 34½" W x 25"D WINDOW: 60" H x 30" W

## **683482** Wide Mobile Barrier

SIZE: 72" H x 52½" W x 25"D WINDOW: 36" H x 48" W



## 683465 **Small Window Mobile Barrier**

SIZE: 721/4" H x 521/2" W x 25"D WINDOW: 12" H x 12" W

## 683468 **Medium Window Mobile Barrier**

SIZE: 74" H x 521/2" W x 25"D WINDOW: 12" H x 24" W

## 683460 Standard Mobile Barrier

SIZE: 72" H x 341/2" W x 25"D WINDOW: 24" H x 30" W

## 683484 Regular Mobile Barrier

SIZE: 641/2" H x 401/2" W x 25"D WINDOW: 24" H x 36" W

## 683458 "Shorty" Mobile Barrier

SIZE: 45" H x 341/2" W x 25"D WINDOW: 24" H x 30" W

## RACKS

## Save Space + Increase Apron Life

INFAB offers apron racks for storing lead aprons used in healthcare facilities. These racks are sturdy and can hold multiple aprons, ensuring easy access for staff. They come in a variety of sizes and designs.






10 Hanger Mobile

Apron Valet

48½ " 683438 20 Hanger Mobile Apron Valet



683419 Mini Mobile Apron Rack with Hangers



683407 Deluxe Mobile Apron Rack



5 Arm Budget Saver Mobile Apron Rack



683427 10 Arm Budget Saver Mobile Apron Rack

## **ACCESSORIES**

## Gloves, Cleaners, Disposables, + More

INFAB offers a variety of options for radiation protection gloves and cleaning solutions to meet your equipment cleaning needs. Our wide range of products allows you to choose the best solution for your specific task.



**SPLIT PALM MITTENS WITH FLAP** 18 Colors Available



REVOLUTION MAXI-FLEX 5 FINGER LEAD GLOVES

COLOR: Black, Red, Blue (shown) and 15 additional colors



## **GLOVES**

INFAB offers customizable gloves in various colors, styles, and lead equivalencies. Choose from reusable or disposable with heavy, medium, or light attenuation level rating.





REVOLUTION **GLOVES** 



**LEAD-FREE NEOPRENE RADIATION PROTECTION GLOVES** 

## **CLEANERS**

INFAB's Scrubbles apron cleaner removes barium, blood, and other radiology stains while leaving behind a clean citrus scent. Additionally, Clorox Healthcare products are also available for cleaning and disinfection in healthcare facilities



Infab is proud to be a global cleaning solutions partner with Clorox.

- · No harsh chemical odor or fumes
- · No fragrance added
- · Nonflammable, unlike quaternary ammonium-alcohol-based wipes
- · No volatile organic compounds (VOCs)









## **CASSETTE COVERS**

Our cassette covers provide exceptional protection against contamination during radiographic process. They protect film, imaging plates and intensify screens.





## Avaiable in a Variety of Sizes











## **CURTAINS**

INFAB's lead curtains are customizable and come in any size, color, or fabric. They have lead equivalency options of .25mm, .50mm, and .75mm.





Straight Track



Curved Track

## **DENTAL**

Our dental aprons protect patients and staff members with colorful designs and the ability to customize with patterns and logos.







DTL-PANO Pano Cape Dental Apron

## **PATIENT**

Our eye, thyroid, ovarian, and gonad protective shields ensure optimal patient care during medical procedures.



PP-ES Disposable Eye Shields



PP-TS Disposable Thyroid Shield



HSV **Half Apron** 



683400 **Patient Protection Set** 



684120 **Ovarian Shields** 



PGS **Gonad Shields** 

## Sterile Drapes

INFAB's sterile drapes are high-quality, sterilized materials for surgical procedures. They come in various sizes and design. They provide a barrier to reduce contamination and protect both the patient and equipment during surgical procedures.





C ARM DRAPE

C ARM DRAPE KIT



FOOT SWITCH COVER

Contact Us M-F 6:30am-5pm PDT info@infabcorp.com www.lnfabCorp.com

Infab Corporation 1040 Avenida Acaso Camarillo, CA, USA 93012 Phone: 805-987-5255 Fax: 805-482-8424







